



FLOW CONTROL SWITCH

Confidently Control Fluid Flow

FLOW CONTROL SWITCH

The Merit Medical Flow Control Switch gives you the confidence you need when stopping and starting flow during critical procedures. This feature rich, integrated, one-handed, single channel device is designed with clinician and patient safety in mind.



FEEL THE DIFFERENCE

An internal closing mechanism allows the device to be locked in the off position without any sharp or protruding parts that could pinch or tear your glove resulting in waste, inconvenience and possible exposure.

HEAR THE DIFFERENCE

Unlike the competition, the Merit Medical Flow Control Switch utilizes a loud audible click when the valve is closed and locked, giving you confidence that fluid flow has been stopped.

EXPERIENCE THE DIFFERENCE

The ergonomically designed, low-profile, non-slip surface gives you more control without rotating an attached device.

Catalog Number	Description	Box Quantity
U1SHC	Flow Control Switch	24

TECHNICAL SPECIFICATIONS

Body Material:	Polycarbonate
Valve:	Silicone
Slider Switch:	Acetal
Pressure Rating:	Open 1200psi Closed 174psi
Lumen I.D.	0.065 inches
Connectors:	Male to Female

Ask today to evaluate the Merit Medical Flow Control Switch.



Merit Medical Systems, Inc. • 1600 West Merit Parkway • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT

Merit Medical Europe • Horsterweg 24 • 6199 AC MAASTRICHT- AIRPORT The Netherlands • +31 43 35 88 222

Merit Medical Ireland Ltd. • Parkmore Business Park West Galway, Ireland • +353 (0) 91 703 733

Free phone from specific country:

Belgium 0800 72906 • France 0800 91 60 30 • Germany 0800 182 0871 • Luxembourg 0800 22522 • Netherlands 0800 022 8184 • UK 0800 973 115

www.merit.com

402244001/A ID 070111