



Surgical Stabilization Products

*Assistant™ Attachment with StableSoft™ Technology
Hercules 3™ Stabilizing Arm*



Simple. Smart. Solutions.

Terumo's unique Surgical Stabilization Products offer positive positioning and retraction during a variety of surgical procedures where a person's hand is traditionally needed.

Surgical Stabilization Products

Our Assistant™ Attachment Brings a Smart Solution to the Table.

Assistant™ Attachment with StableSoft™ Technology

Available in 4 sizes — 2 and 3 fingers, long and short — to accommodate various anatomy and surgical applications. Because they bend to conform to anatomical shapes, they are ideal for procedures requiring a human hand. The textured, silicone-coated surface provides exceptional grip and positional stability, helping to maximize field of view.



2L-2 Finger Long

1 ¼" (3.2 cm) W x 3 ½" (8.9 cm) L
Shaft length: 1 ½" (3.8 cm)



2S-2 Finger Short

1 ¼" (3.2 cm) W x 2" (5.1 cm) L
Shaft length: 1 ½" (3.8 cm)



3L-3 Finger Long

2 ⅛" (5.4 cm) W x 3 ½" (8.9 cm) L
Shaft length: 1 ½" (3.8 cm)



3S-3 Finger Short

2 ¼" (5.7 cm) W x 2" (5.1 cm) L
Shaft length: 1 ½" (3.8 cm)

Hercules™ Stabilizing Arm

Reusable stainless steel Hercules Stabilizing Arms reduce supply costs and environmental impact. Universal attachments facilitate ease of use on most sternal retractors while the Arm's low profile articulation maximizes flexible positioning and strength.



Ordering Information

Assistant Attachment with StableSoft Technology

CATALOG #	DESCRIPTION	UNITS/CASE
001-401-209S	Assistant Attachment with StableSoft Technology, 3L Single use, malleable, three fingers, long	5
001-401-210S	Assistant Attachment with StableSoft Technology, 3S Single use, malleable, three fingers, short	5
001-401-211S	Assistant Attachment with StableSoft Technology, 2L Single use, malleable, two fingers, long	5
001-401-212S	Assistant Attachment with StableSoft Technology, 2S Single use, malleable, two fingers, short	5

Hercules Stabilizing Arm

CATALOG #	DESCRIPTION	UNITS/CASE
001-401-161	Stabilizing Arm, Hercules 3, Universal Top mount, low profile, reusable stainless steel	1



For more information on this or other Terumo cardiovascular products, contact:

TERUMO CARDIOVASCULAR GROUP
6200 Jackson Road
Ann Arbor, MI 48103-9300
USA
734.663.4145
734.663.7981 fax
www.terumo-cvgroup.com

TERUMO EUROPE N.V.
Authorized EC Representative
Interleuvenlaan 40
3001 Leuven
Belgium
32.16.38.12.11
32.16.40.02.49 fax
www.terumo-europe.com

TERUMO EUROPE N.V.
CARDIOVASCULAR DIVISION
Ludwig-Erhard-Strasse 6
65760 Eschborn
Germany
49.6196.80.230
49.6196.8023.555 fax
www.terumo-europe.com

TERUMO LATIN AMERICA CORPORATION
8750 NW 36th Street, Suite 600
Miami, FL 33178
USA
305.477.4822
305.477.4872 fax

TERUMO CORPORATION
44-1, 2-chome
Hatagaya, Shibuya-Ku
Tokyo, 151-0072
Japan
81.3.3374.8111
81.3.3317.6010 fax
www.terumo.co.jp