



Specifications

Product Code	Catheter Size	Gauge Size	Length(cm)							Spring Guidewire (J-Tip)			Dilator	Introducer Needle	
			6	8	10	13	16	20	30	Size				18Gx2 1/2"	20Gx1 1/2"
										0.018"	0.025"	0.035"			
SINGLE LUMEN															
CV-401-16	4F	18G					✓				✓		5F		✓
CV-401-20	4F	18G						✓			✓		5F		✓
CV-501-16	5F	16G					✓				✓		6F	✓	
CV-501-20	5F	16G						✓			✓		6F	✓	
CV-701-16	7F	14G					✓				✓		7F	✓	
CV-701-20	7F	14G						✓			✓		7F	✓	
DOUBLE LUMEN															
CV-402-06	4F	22/22G	✓								✓		5F		✓
CV-402-10	4F	22/22G		✓							✓		5F		✓
CV-402-16	4F	22/22G					✓				✓		5F		✓
CV-502-10	5F	18/18G		✓							✓		6F		✓
CV-502-16	5F	18/18G			✓						✓		6F		✓
CV-502-20	5F	18/18G					✓				✓		6F		✓
CV-702-16	7F	16/16G					✓				✓		8F	✓	
CV-702-20	7F	16/16G						✓			✓		8F	✓	
CV-702-30	7F	16/16G							✓		✓		8F	✓	
CV-802-16	8.5F	13/13G					✓				✓		10F	✓	
CV-802-20	8.5F	13/13G						✓			✓		10F	✓	
CV-802-30	8.5F	13/13G							✓		✓		10F	✓	
TRIPLE LUMEN															
CV-503-08	5.5F	20/22/22G		✓							✓		6F		✓
CV-503-13	5.5F	20/22/22G			✓						✓		6F		✓
CV-503-30	5.5F	20/22/22G							✓	✓			6F		✓
CV-703-10	7F	16/18/18G		✓							✓		8F	✓	
CV-703-16	7F	16/18/18G					✓				✓		8F	✓	
CV-703-20	7F	16/18/18G						✓			✓		8F	✓	
CV-703-30	7F	16/18/18G							✓		✓		8F	✓	
CV-713-10	7F	14/18/18G		✓							✓		8F	✓	
CV-713-16	7F	14/18/18G					✓				✓		8F	✓	
CV-713-20	7F	14/18/18G						✓			✓		8F	✓	
CV-713-30	7F	14/18/18G							✓		✓		8F	✓	
CV-803-16	8.5F	14/18/18G					✓				✓		10F	✓	
CV-803-20	8.5F	14/18/18G						✓			✓		10F	✓	
CV-803-30	8.5F	14/18/18G							✓		✓		10F	✓	
QUAD LUMEN															
CV-804-16	8.5F	14/16/18/18G					✓				✓		10F	✓	
CV-804-20	8.5F	14/16/18/18G						✓			✓		10F	✓	
CV-804-30	8.5F	14/16/18/18G							✓		✓		10F	✓	


Optional Features

Family/Configuration	Description of Configuration	Remark
CV-XXX-XX	CVC with Spring Guidewire and straight needle	NA
CVN-XXX-XX	CVC with Kink-Resistance Guidewire and straight needle	N-Refer to Kink-Resistance Guidewire
CVNR-XXX-XX	CVC with Kink-Resistance Guidewire and Raulerson Syringe	R-Refer to Raulerson Syringe
CV-XXX-XX Y	CVC with Spring Guidewire and Y-Valve Needle	Y-Refer to Y-Valve Needle
CV-XXX-XX T	CVC kit with Procedure Tray	T-Refer to Procedure Tray
CVR-XXX-XX	CVC with Raulerson Syringe	R-Refer to Raulerson Syringe

All complete CVC kit includes Extension Line Clamp, Injection Site, 1 ea. of Luer Slip Syringe, 1 ea. of #11 Short Scalpel, 1 ea. of Suture Wing with Clamp Fastener & 1 ea. of Needle Sponge.

CENTRAL VENOUS CATHETER

Multi-lumen Catheter



Bioptimal™

BIOPTIMAL INTERNATIONAL PTE. LTD.

36 Jalan Tukang
Singapore 619266
Tel: +65 6213 5777
Fax: +65 6213 5737
Email: sales@bioptimalg.com



www.bioptimalg.com

CENTRAL VENOUS CATHETER

Multi-lumen Catheter

Always convenient.
Always ready.



Soft Tip

Superior geometric soft tip design makes insertion safe and easy.

Polyurethane Material

Special formulated biocompatible polyurethane material, provides strength during insertion, yet softens in body.

Catheter Marking

Radiopaque material with clear, definite marking, ensures correct placement of the catheter.



Suture Wing & Clamp Fastener

Keep catheter in place thereby reducing the risk of mechanical phlebitis and irritation.

Quality Components for Optimal Performance



Tegaderm Transparent Dressing (optional)

Transparent window provides visibility, stabilizes catheter and protects from outside contamination.

Introducer Needle

Available in Straight or Y needle. Y valve needle designed to prevent back-flow of blood and air embolism.

Advancer

Single hand operated advancer helps to straighten the "J" Tip spring guidewire and advancing it to the proper position.

Guidewire

Specially designed guidewire with a soft flexible tip to prevent vessel perforation (vascular trauma), yet provides good torque to ensure firm insertion. Markings on guidewire provide a good insertion indication.

Needle Sponge

Specialized sponge for holding used needle preventing the potential of needle piercing.

Contents of the CVC Procedure Tray

- Needle 22G x 1 1/2"
- Needle 25G x 5/8"
- Paper Towel 17" x 22"
- Gauze White
- Fenestrated Drape
- 5cc Syringe with Luer Slip
- Sponge Swabs
- Suture: Braided Silk 3-0 with Straight Needle
- CVC Tray-Medication Cup
- CVC Kit
- Blue Wrap (70 x 70)
- Tegaderm Transparent Dressing (optional)

The CVC Procedure Tray is optional and can be purchased together with the standard central venous catheter kit as a complete set. Please approach our sales representatives for further enquires.

Biopital's Basic CVC Kit

- Polyurethane Catheter
- "J" Tip Guidewire complete with Dispenser & Advancer
- Introducer Needle (Straight Needle or Y Valve Needle)
- Tegaderm Transparent Dressing (optional)
- Suture Wing & Clamp Fastener
- Extension Line Clamp
- Needle Sponge
- Vessel Dilator
- 5cc Syringe with Luer Slip
- Scalpel #11, Short
- Injection Site

Customization is available, please contact your local sales representative for more information.

